## ICP – Mass Spectrometry

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Yongbo Dan<sup>1,2</sup>, Weilan Zhang<sup>3</sup>, Xingmao Ma<sup>2,3</sup>,  
Honglan Shi<sup>1,2</sup>, Chady Stephan<sup>4</sup>

<sup>1</sup>Department of Chemistry, Missouri University of  
Science and Technology

<sup>2</sup>Center for Single Nanoparticle, Single Cell, and  
Single Molecule Monitoring (CS<sup>3</sup>M), Missouri  
University of Science and Technology

<sup>3</sup>Department of Civil and Environmental  
Engineering, Southern Illinois University

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Table 1. NexION 300/350D Instrumental and Analytical Parameters.

Parameter	Value
Nebulizer	Concentric (glass)
Nebulizer Flow	1.08 L/min
Spray Chamber	Baffled Cyclonic (glass)
ICP RF Power	1600 W
Analyte	Au
Mass	197 amu
Dwell Time	0.1 ms
Settling Time	0 ms
Sampling Time	100 sec
Number of Data Points Acquired	1 million per sample
Au Density	19.3 g/cm <sup>3</sup>

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
Figure 1. Particle size distribution histogram of enzyme-treated 50 nm AuNP (without plant tissue).

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
Figure 2. (a) Raw data for reagent blank (reagent blank: enzyme in 2 mM citrate solution, without plant tissues and Au NPs); (b) Raw data for control plant without exposure to Au NPs; (c) Raw data for spiking  $4.7 \times 10^4$  NPs/mL of 100 nm Au NPs into control plant sample; (d) Size distribution histogram for spiking  $4.7 \times 10^4$  /mL of 100 nm Au NPs into control plant sample.

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
Figure 3. (a)&(b) Raw data of duplicated tomato plants exposed to 5 mg/L of 40 nm Au NPs for 4 days; (c) Size distributions histograms of duplicated tomato plants exposed to 5 mg/L of 40 nm Au NPs from Figure 4(a) and 4(b); (d) Size distribution histogram of spiking  $4.7 \times 10^4$  particles/mL of 100 nm Au NPs into tomato plants exposed to 5 mg/L of 40 nm Au NPs.

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1. Yongbo Dan, Weilan Zhang, Runmiao Xue, Xingmao Ma, Chady Stephan, Honglan Shi, 2015, "Characterization of Gold Nanoparticles Uptake by Tomato Plants Using Enzymatic Extraction Followed by Single Particle Inductively Coupled Plasma-Mass Spectrometry Analysis", *Environmental Science and Technology*, 49(5):3007-3014.

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Component	Part Number
Green/Orange Flared Peristaltic Pump Tubing (package of 12)	N0777042
Gold Nanoparticles in Water (30 nm)	N8142300
Gold Nanoparticles in Water (50 nm)	N8142302
Gold Nanoparticles in Water (80 nm)	N8142305
Gold Nanoparticles in Water (100 nm)	N8142307
Autosampler Tubes	B0193233 (15 mL) B0193234 (50 mL)

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版权所有 ©2014, PerkinElmer, Inc. 保留所有权利。PerkinElmer® 是PerkinElmer, Inc. 的注册商标。其它所有商标均为其各自持有者或所有者的财产。